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Instructions for Use
Istruzioni per l'uso
Manuel d'utilisation

Sono*One*

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This instruction for use is an integral part of the device. It must be stored with the device and kept accessible at all times for anyone who is permitted to operate the device. The instructions for use are valid as of June 2022.

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SonoOne is an ultrasound therapy device for application of ultrasound in order to heat tissue for purposes of treatment of diseases of the musculoskeletal system.

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SonoOne is not intended for combination with other devices.

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Heating tissue in cases of diseases of the musculoskeletal system e.g. Epicondylopathy, Joint diseases.

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EN

- Vague pain symptoms
- Diseases where heat should not be used, for example, acute inflammatory diseases
- Diseases where mechanical impacts are contraindicated, e.g. deep vein thrombosis
- Restricted blood flow
- Suspected diseases of the cardiovascular system
- Haemorrhagic diathesis
- Cervical vertebrae – do not use ultrasound higher than C 3
- Ultrasound treatment of parenchymatous or heat-sensitive organs (testes, eyes, gravid uterus, liver, kidneys etc.)
- Anaesthetised skin areas
- Temperature sensitivity disorders
- After treatment with ionizing radiation
- Epiphyseal plates
- Tumours
- Electronic pacemakers
- Application on or in the vicinity of open wounds
- Impaired reflexes or pain sensitivity
- Application to the abdomen, pelvis, or lumbar area when pregnant or presumed pregnant

Metal implants and endoprotheses: Dynamic ultrasound radiation in low dosage is nowadays no longer a matter for concern. No side effects are known in the therapeutic application of ultrasound.

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To avoid the risk of electric shock, disconnect the device from the mains supply by pulling the plug before carrying out any maintenance or cleaning work.

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Before use, ensure that the device is powered via a properly earthed plug with a grounded outlet (electrical installation according to DIN VDE 0100 Part 710). The device must only be operated with the supplied power cord. The power cord must be protected against mechanical stress.

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Magnetic and electric fields can influence the function of the device. Therefore, the device must not be operated in the vicinity of devices that generate strong electromagnetic fields (X-ray or diathermy devices, MRI scanners). Please keep a safe distance of several metres.

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The device is not suitable for use in areas with an explosive, combustible or combustion-promoting atmosphere.

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During use, the device must be placed in such a way that direct access to the device's central power supply is possible, so that it can be disconnected from the mains at any time.

Inspect the device and the ultrasound head before use. If damaged, it must not be used. Do not allow sharp or pointed objects to come into contact with the ultrasound transducer since the aluminium head can be easily scratched.

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The use of coupling agents other than the special SonoPlus ultrasound gel can damage the ultrasound transducer.

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Only accessories provided by Zimmer MedizinSysteme GmbH are to be used.

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The device may cause malfunctions or may interfere with the operation of equipment in the vicinity. It may be necessary to take appropriate remedial action, such as realignment, re-arrangement of the device or installation shielding.

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During the service life of the device, no changes may be made to the device or medical system.

The entire medical system is suitable for use in the vicinity of the patient.

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In the event of regular use of the ultrasound head in the water bath, the user must take protection measures. All body parts of the user immersed in the water must be protected by appropriate measures. Suitable protective measures are so-called "aerated clothes".

This can be neoprene gloves or latex over cotton gloves.

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This device is intended to be used exclusively by medical professionals.

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The use of this product with settings or for applications other than those specified in the user manual may lead to hazards resulting from the uncontrolled release of ultrasound energy.

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The patient must not be left unattended during therapy.

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Any treatment instructions regarding treatment location, duration and intensity require medical knowledge and should be given by authorised physicians, therapists and medical paraprofessionals. These instructions must be followed.

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Use in wet areas is not permitted and may in case of non-compliance lead to considerable damage and endanger both the patient and the user.

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Dispose of the packaging material properly. Make sure that it is not accessible to children.

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In cases of patients with implants or an implanted electronic device do not conduct treatment before the risk has been assessed and found negligible.

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Patients may not be connected to a high-frequency surgical device at the same time. This can lead to burns.

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Fig. 4 Parameter screen

Fig. 4 Videata Parametri

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Fig. 5 Switches and connector sockets

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Description of screens, displays, windows, buttons

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Note: Activate all windows and buttons by pressing directly on the screen with your finger.

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Mounting the cables / therapy instruction in brief

Connect the mains cable with the connection for the cable. Connect the ultrasound head with the connecting socket for the ultrasound head. Switch on the device. After the program screen appears choose your programme, adjust intensity and the therapy time and press "Start". Therapy ends when therapy time has elapsed.

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EN	Mains voltage	100 - 240 V / 50 Hz / 60 Hz 220 V / 60 Hz
	Power Consumption	max. 60 VA
	Protection class	Class II
	Applied Part	Type BF, Transducer ultrasound head
	Operation Mode	Continuous operation
	Dimensions	238 (L) x 138 (W) x 95 (H) mm
	Weight	1kg
	Operation	10°C to +40°C, 30% to 75% relative humidity without condensation, at 800 hPa - 1060 hPa
	Transport/Storage	-10°C to +50°C, 10% to 90% relative humidity without condensation, at 700 hPa - 1060 hPa
	Ultrasound head	5 cm ² , A _{ER} = 3.0 cm ² BN ≤ 5.0 IPX 7: protected from immersion in water with a depth of up to 1 m for up to 30 min.
	Frequency	1 MHz and 3MHz
	Intensity levels	0.1 to 3 W/cm ² at intervals of 0.1 W/cm ²
	Maximum Power	6 W _{eff} in continuous mode, 9 W _{peak} in pulse mode
	Accuracy	± 20%
	Ultrasound types	1. Continuous ultrasound 2. Pulsed ultrasound, adjustable pulse frequencies: 10Hz, 20 Hz, 50 Hz, 100 Hz Duty ratio: 1: 1, 1: 2, 1: 5, 1: 10
	Interchangeability	Ultrasound heads are calibrated ex works and can be easily replaced. Storage and transport in original packaging only.

Subject to technical modifications!

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Make sure that during cleaning or disinfection no liquid penetrates the device. Do not use sprays.

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If during cleaning or disinfecting liquid penetrates the device, please put the device out of service, protect it from being used again and contact your service representative.

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Make sure that when cleaning and disinfecting the labelling of the device (such as warnings, labels, identification plate) is not damaged.

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The device and its applied part are considered as uncritical in relation to hygiene due to the use on non-injured and healthy skin.

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Cleaning: In the event of visible contamination, the housing and all cables can be cleaned using commercially available alcohol-free plastic cleaners. Wipe the surface until the dirt is removed, using a soft cloth soaked according to the specifications of the manufacturer of the cleaning agent but not dripping wet.

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Disinfection: We recommend that disinfection is to be carried out at least once a week, as well as if there is any indication of contamination. Consult with your hygiene specialist when doing so. Always perform cleaning prior to disinfection.

Housing and all cables can be disinfected using disinfectant wipes. Use a commercially available alcohol-free disinfectant for metal and plastic, with bactericidal, virucidal and fungicidal properties. Observe the instructions for use of the respective manufacturer. Wipe all surfaces using a soft cloth soaked according to the specifications of the manufacturer of the disinfectant, but not dripping, or with a cloth pre-impregnated with disinfectant (wipes). Also observe requirements for drying or post-cleaning, where applicable.

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If flammable solutions are used for cleaning and disinfecting, sufficient time must be allowed for the solution to evaporate before using the device. Otherwise, it may constitute a fire hazard.

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Use the device only in a hygienic environment.

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EN The devices have a CE mark

in accordance with the EC directive on medical devices 93/42/EEC.

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EN Error Code list

224/1	Communication Error		
402/2	Ultrasound applicator has overheated	407/2	PLL fault
406/2	Ultrasound head not connected	410/2	Output Voltage abnormal

In some cases switching the device on an off may solve the problem. Service and repair only by trained staff. Service manual available.

Functional test

Cover ultrasound head with coupling gel, set intensity, start therapy. Coupling display must indicate over 90 %.

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EN For questions or in case of equipment malfunctions, please contact the address listed:

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EN Sono *One* is not listed in Annex 1 of MPBetreibV. The device is not listed in Annex 2 of MPBetreibV. In Germany, the DGUV regulation 3 (Electrical Systems and Equipment) in its current version must also be observed.

For disposal, the device may only be returned to Zimmer MedizinSysteme in its original packaging. In foreign countries (outside of Germany) please refer to national regulations for disposal. Contact your distributor if necessary.

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These notes apply to the operation of the device in Germany. If applicable, observe any divergent national regulations that apply in your country.

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Medical electrical equipment such as the SonoOne are subject to special precautions with regard to EMC (electromagnetic compatibility) and must be installed and commissioned in accordance with the EMC information contained in the user manual or the accompanying documents.

Portable and mobile RF communications equipment (e.g. mobile phones) can affect electrical medical equipment.

SonoOne may only be operated with the original parts specified in the list of parts and accessories included within the scope of delivery. The operation of the device with other parts may result in increased emissions or in the device's decreased interference immunity.

Guidelines and Manufacturer's Declaration - Electromagnetic Emissions		
The SonoOne device is intended for use in the electromagnetic environment specified below. The customer or the user of the SonoOne must ensure that it is operated in such an environment.		
Emissions measurements	Compliance	Electromagnetic environment - guidelines
RF emissions in accordance with CISPR 11	Group 1	The SonoOne device uses RF energy solely for its internal functioning. Its RF emission is therefore very low and it is unlikely that this will cause interference to neighbouring electronic Rehab-series
RF emissions in accordance with CISPR 11	Class A	The SonoOne device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Emissions of harmonics in accordance with IEC 61000-3-2	Category A	
Emissions of voltage fluctuations / flickers in accordance with IEC 61000-3-3	Compliant	

Table 201 according to EN 60601-1-2: 2006-10

The device must not be used in direct proximity to or stacked directly on top of another similar device. If operation near to or stacked on top of another device is unavoidable, the device should be monitored to verify its proper operation within this setup.

Guidelines and manufacturer's declaration - electromagnetic immunity			
The SonoOne device is intended for use in an electromagnetic environment, as specified below. The customer or the user of the SonoOne device should ensure that it is used in such an environment.			
Immunity tests	IEC 60601 - Test Level	Compliance level	Electromagnetic environment - Guidelines
Electrostatic discharge (ESD) in accordance with IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be made from wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient / burst in accordance with IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output cables	± 2 kV for mains cables Not applicable	The supply voltage quality must correspond to that of a typical commercial or hospital environment.
Surges in accordance with IEC 6100-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Not applicable The supply voltage quality must correspond to that of a typical commercial or hospital environment.
Voltage dips, short-term interruptions and voltage variations in accordance with	< 5% UT (> 95% dip in UT for ½ period) 40% UT (60% dip in UT for 5 periods) 70% UT (30% dip in UT for 25 periods) < 5% UT (> 95% dip in UT for 5 seconds)	< 5% UT (> 95% dip in UT for ½ period) 40% UT (60% dip in UT for 5 periods) 70% UT (30% dip in UT for 25 periods) < 5% UT (> 95% dip in UT for 5 seconds)	The supply voltage quality must correspond to that of a typical commercial or hospital environment. If the user of the SonoOne requires continued operation, even in the case of interruptions in the power supply, it is recommended that the SonoOne be powered from an uninterruptible power supply or a battery.
Magnetic field of the supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at mains frequency should have the typical values found in a business or hospital environment.
Note: UT is the mains AC voltage before application of the test level			

Key features of the SonoOne are: interference-free delivery ultrasound, with appropriately configured parameters and trouble-free operation of all functions.

Guidelines and manufacturer’s declaration - electromagnetic immunity

The SonoOne device is intended for use in an electromagnetic environment, as specified below. The customer or the user of the SonoOne device should ensure that it is used in such an environment.

Immunity tests	IEC 60601 - Test Level	Compliance level	Electromagnetic environment - Guidelines
Conducted RF disturbance variables in accordance with IEC 61000-4-6	3 V _{RMS} 150 kHz to 80 MHz	3 V _{RMS} 150 kHz to 80 MHz	Portable and mobile RF communications equipment should not be used any closer to the SonoOne including cables, than the recommended separation distance calculated using the equation relevant for the transmission frequency. Recommended separation distance: $d = 1.17 \sqrt{P}$ $d = 0.35 \sqrt{P}$ for 80MHz to 800 MHz $d = 0.7 \sqrt{P}$ for 800 MHz to 2.5 GHz where P is the power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). In accordance with an on-site investigation, for all frequencies, the field strength of a stationary radio transmitter should be less than the coincidence level. Errors may occur in the vicinity of devices marked with the following symbol: 
Radiated RF disturbance variables in accordance with IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m 80 MHz to 2.5 GHz	

NOTE 1: At 80 Hz and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) The field strengths of fixed transmitters e.g. base stations of mobile telephones and mobile radio devices, amateur radio, AM and FM radio broadcasting, and TV broadcasting cannot be accurately predicted. In order to assess the electromagnetic environment with regard to fixed RF transmitters, a study of the electromagnetic phenomena at the location should be considered. If the measured field strength in the location in which the SonoOne device is used exceeds the above compliance level, the SonoOne device must be monitored to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as the realignment or relocation of the SonoOne device.
- b) In the frequency range of 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the SonoOne device

The SonoOne device is intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or the user of the SonoOne device can thus help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SonoOne device – depending on the output power of the communications equipment as indicated below.

Rated power of the transmitter W	Separation distance depending on the frequency of transmitter m		
	150 KHz to 80 MHz $d = 1.17 \sqrt{P}$	80–800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.12	0.04	0.07
0,1	0.37	0.11	0,22
1	1,17	0,35	0,7
10	3,70	1,11	2,21
100	11,67	3,5	7,0

For transmitters for which the maximum nominal output is not specified in the above table, the recommended separation distance d in metres (m) can be determined using the equation given for the corresponding column, where P is the maximum rated power of the transmitter in watts (W) in accordance with the manufacturer of the transmitter.
 NOTE 1: At 80 and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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